

Lapiplasty®

Sterile-Packed Release Instruments

SpeedRelease™

Guided Release Instrument

Sterile-packed, single-use instrument designed for quick and controlled release of the sesamoidal suspensory ligament and other soft tissues.

- **Guided tip** to direct insertion within the lateral joint capsule
- **Cutting edge** for quick and controlled release of the contracted soft tissue
- **Sterile-packed** for convenient delivery and consistent sharpness



TriTome™

Triple-Edge Release Instrument

Sterile-packed, single-use instrument designed to release between the metatarsal bases for the Adductoplasty™ Procedure and other applications.

- **Three cutting edges** for quick and controlled soft-tissue release
- **Thin 1.5mm cutting end** to access challenging anatomy
- **Sterile-packed** for convenient delivery and consistent sharpness




TREACE
Medical Concepts, Inc.

The Leader in Hallux Valgus Surgery™


1st MTP Lateral Release with the SpeedRelease™ Instrument

Key Surgical Steps




1. Lateral Capsule Incision

Make a small vertical incision in the lateral capsule of the 1st MTP joint.



2. Creation of Capsular Pocket

Insert a hemostat into the lateral capsular incision to create a soft-tissue pocket.



3. Insertion of SpeedRelease™ Instrument

Insert the blunt tip of the SpeedRelease™ instrument into the pocket in the lateral capsule.



4. Sesamoidal Ligament Release

Advance plantar and posteriorly, between the sesamoids and metatarsal head, to release the sesamoidal suspensory ligament.

Ordering Information

SN20 SpeedRelease™ Guided Release Instrument

SN21 TriTome™ Triple-Edge Release Instrument

Refer to Treace Medical Concepts LBL 09-00001L for full Instructions for Use.



To learn more about benefit and risks, visit Lapiplasty.com

The Leader in Hallux Valgus Surgery™

TREACE
Medical Concepts, Inc.