

# FastGrafter<sup>®</sup>

Autograft Harvesting System



**Fully-sterile system for quick and efficient  
harvest of autogenous bone**

**TREACE<sup>®</sup>**  
Medical Concepts, Inc.

*The Leader in Hallux Valgus Surgery™*

# FastGrafter®

## Autograft Harvesting System


The **FastGrafter® Autograft Harvesting System** is a sterile-packed, single-use device designed for quick and efficient harvest of cancellous autogenous bone from the calcaneus, distal tibia, and other harvest sites through a minimal incision approach.

- **Single-piece** harvester designed to reduce instrumentation and system complexity
- **Morselizing cutting** tip penetrates cortex and morselizes bone during harvest
- **Sterile-packed system** designed for quick and efficient harvest of autograft bone



# Calcaneal Autograft Harvest

## Key surgical steps\*



1.

### Incision and Dissection


Make a small incision over the lateral aspect of the calcaneus, posterior and inferior to the peroneal tendon and sural nerve, approximately 3cm below and 2cm posterior to the distal fibula. Use blunt dissection to expose the bone.



2.

### Bone Graft Harvest

Insert the **FastGrafter® Harvester** into AO attachment on a powered driver. Place the cutting tip of the **Harvester** onto the exposed bone surface. Beginning at low speed in the forward direction, advance the **Harvester** to the desired depth.



3.

### Additional Bone Graft Harvest

Reinsert the **Harvester** tip into the original harvest site and make a pass approximately 20-30 degrees from the original harvest path.




4.

### Removal of Morselized Bone Graft

Disassemble the **Harvester** from the powered driver. Over a sterile container, insert the **Pusher** through the distal tip of the **Harvester** to expel the morselized graft through the proximal opening of the **Harvester**.

# Additional Harvest Sites

The **FastGrafter® Autograft Harvesting System** can be used percutaneously through multiple harvest sites including the calcaneus, proximal tibia, and distal tibia.



## Ordering Information

**SK27** FastGrafter® Autograft Harvesting System (7mm)

Learn more about the **FastGrafter® Autograft Harvesting System** at [Laplasty.com](https://Laplasty.com)



*The Leader in Hallux Valgus Surgery™*

**TREACE**  
Medical Concepts, Inc.

The content herein is for informational purposes and intended for healthcare professionals only. Treace Medical Concepts, Inc. does not dispense medical advice. A surgeon should rely on his or her own professional clinical judgment when deciding whether to use a specific product when treating a particular patient. Prior to use of the device, the surgeon should refer to the instructions for use for complete indications, contraindications, warnings, and risks. See FastGrafter IFU (LBL 1507-9005). For complete risk and benefit information, visit [Treace.com](https://Treace.com).