


Percuplasty™

Percutaneous 3D Bunion Correction®



Control and Efficiency at Every Step

TREACE
Medical Concepts, Inc.
A Step Ahead.™

Percuplasty™

Percutaneous 3D Bunion Correction®

- **Precise 3D Correction and Targeting**


Elegantly instrumented to dial in and maintain hands-free triplanar correction and accurately target wire placement

- **Throw Wires Without Skiving**

Robust 1.4/1.6mm cobalt chrome wires with a laser sharp tip designed to reduce skiving

- **Rapid Screw Insertion**

Self-drilling beveled screws help simplify the workflow and save valuable time in the OR



Ordering Information

3.0mm Screws		4.0mm Screws	
SF30XX (XX = Length) Available Lengths: 18mm – 50mm In 2mm increments		SF40XX (XX = Length) Available Lengths: 24mm – 64mm In 2mm increments	
Shannon Burr	Shannon Burr	Wedge Burr	
TM12013 2mm x 13mm	TM12220 2.2mm x 20mm	TM22913 3mm	

Before use of the system, the surgeon should refer to the appropriate instructions for use and surgical technique for complete warnings, precautions, indications, contraindications, and adverse events. Risks include, but are not limited to: infection, pain, discomfort from the presence of the implant, loosening of the implant, and loss of correction with nonunion or malunion. If any of these occur, additional treatments may be needed. Additional information about risks, warnings, and instructions is available at Lapiplasty.com/surgeons/labeling.



**To learn more,
visit Treace.com**